



Be Right™



Electrode Cleaning Solution for Proteins/Organics Samples, 250 mL

Product #:

C20C370

ZAR Price:

Contact Hach

Ships within 1 week

 **Hazardous**

Items with this mark may be considered hazardous under some shipping conditions.

If necessary, we will change your selected shipping method to accomodate these items.

Electrode cleaning solution to remove proteins residues organic contaminations/deposits such as food stuff from pH sensing element. 250 mL Bottle

Specifications

Description:	Remove proteins residues organic deposits
Formula:	PEPSIN in HCl
Quantity:	250 mL
Solution Type:	Cleaning