



Cell Assemblies, Flow-Through

Product #:

4521500

 **Hazardous**

OBSOLETE ITEM

This item is no longer available.

Items with this mark may be considered hazardous under some shipping conditions.

If necessary, we will change your selected shipping method to accommodate these items.

Obtain colorimetric measurements with Hach spectrophotometers and turbidimeters without individual sample cells. Eliminate errors due to scratched or smudged sample cells. No need for optically matched sample cells.

Specifications

Includes:	Assembly tubings, funnels, pour-thru cell
Rinsing volume:	Minimum 20 mL