



Digital titrator cartridge, sulfuric acid, 8.00 N

Product #:

1439101

ZAR Price:

Contact Hach

Ships within 3 days

 **Hazardous**

Items with this mark may be considered hazardous under some shipping conditions.

If necessary, we will change your selected shipping method to accomodate these items.

For base concentration determination by Digital Titrator titration. Digits x 0.01 = meq

Accurate results - when and where you need them

Digital Titrator methods use a high-precision dispensation device and cartridges of concentrated titrant to allow analyses in the lab, the plant or the field.

Fast and Flexible

Interchangeable cartridges allow you to change titrants in seconds. No cleaning or rinsing burets.

Versatile

More than 50 parameter/range combinations are available

Economical

Concentrated titrants allow 40-100 typical titrations per cartridge.

Specifications

Method:	Phenolphthalein End Point
Method Number 3:	10222
Number of tests:	40 - 100
Parameter:	Aqua Ammonia
Platform:	Digital Titrator
Range:	10 - 4000 meq/L Base
Range 2:	5 - 35 % as NH ₃
Titrant:	8.00 N H ₂ SO ₄

Required Accessories

- Gloves, Disposable, Powder Free, Nitrile, Large (Item 2550503)