



# Digital titrator cartridge, hydrochloric acid, 8.00 N

ZAR Price: Contact Hach Ships within 3 days

Items with this mark may be considered hazardous under some shipping conditions.

If necessary, we will change your selected shipping method to accomodate these items.

For base concentration determination by Digital Titrator titration. Digits x 0.01 = meq

#### Accurate results - when and where you need them

Digital Titrator methods use a high-precision dispensation device and cartridges of concentrated titrant to allow analyses in the lab, the plant or the field.

#### Fast and Flexible

Interchangeable cartridges allow you to change titrants in seconds. No cleaning or rinsing burets.

#### Versatile

More than 50 parameter/range combinations are available

#### **Economical**

Concentrated titrants allow 40-100 typical titrations per cartridge.

## **Specifications**

Method: Phenolphthalein End Point

Method Number 2: 8233 (Base)

Number of tests: 40 - 100

Parameter: Acid-Base

Platform: Digital Titrator

Range: 10 - 4000 meq/L Base

Titrant: 8.00 N HCI

### Required Accessories

• Gloves, Disposable, Powder Free, Nitrile, Large (Item 2550503)