



Be Right™



Amino acid reagent powder pillows for LR silica, 25 mL, pk/1000

Product #: 80428

ZAR Price: Contact Hach
No availability shown

For determination of low range silica by the Heteropoly Blue Method (0 - 2.0 mg/L SiO₂). Optional reagent for high range Phosphate by the Amino Acid method. Sample Size: 25 mL. Pack of 1000 Powder Pillows

Specifications

Platform: Powder Pillows

Required Accessories

- Gloves, Disposable, Powder Free, Nitrile, Medium (Item 2550502)
- Gloves, Disposable, Powder Free, Nitrile, X-Large (Item 2550504)
- Gloves, Disposable, Powder Free, Nitrile, Small (Item 2550501)
- Gloves, Disposable, Powder Free, Nitrile, Large (Item 2550503)