



Amino acid reagent powder pillows for LR silica, 25 mL, pk/1000

Product #:

80428

OBSOLETE ITEM This item is no longer available.

Suggested replacements: • 80499

For determination of low range silica by the Heteropoly Blue Method (0 - 2.0 mg/L SiO₂). Optional reagent for high range Phosphate by the Amino Acid method. Sample Size: 25 mL. Pack of 1000 Powder Pillows

Specifications

Platform:

Powder Pillows

Required Accessories

• Gloves, Disposable, Powder Free, Nitrile, Large (Item 2550503)