



Reagent solution, amino acid F, 20 ampules (ultra-low range)

Product #: 2386420

ZAR Price: Contact Hach

Ships within 3 days

For determination of Ultra Low Range silica by the Heteropoly Blue method. Hach method 8282 Range: 3 - $1000 \mu g/L$ SiO2. Ampuled for increased stability

Specifications

Method: Heteropoly Blue

Number of tests: 20

Platform: Solution

Required Accessories

• Gloves, Disposable, Powder Free, Nitrile, Large (Item 2550503)