



Reagent set, aluminum, 50 mL sample, ECR method

Product #:

2443400

ZAR Price:
Available

Contact Hach

 **Hazardous**

Items with this mark may be considered hazardous under some shipping conditions.

If necessary, we will change your selected shipping method to accommodate these items.

For the determination of Aluminium, Eriochrome Cyanine R Method. For 25ml sample.

For determination of Aluminum by the Eriochrome Cyanine R (ECR) Method. Hach Method 8326. Sample Size: 50 mL Range: 0.006 - 0.220 mg/L. Reagent set contains ECR Reagent Powder Pillows, Hexamethylenetetramine Buffer Reagent Powder Pillows and ECR Masking Reagent Solution for approximately 100 tests.

Specifications

Instrument:	DR/3000, DR/2000
Method:	Eriochrome Cyanine R
Number of tests:	100
Range:	0 - 0.220 mg/L Al

Required Accessories

- Gloves, Disposable, Powder Free, Nitrile, Large (Item 2550503)