

Editoria Bases and Base

Calcium Separation Set (for TNT852)

Product #: ZAR Price:

Ships within 1 week

TNT892

Contact Hach

The Calcium Separation Set is required when conducting TNTplus tests for Cadmium (TNT852) in sample matrices that contain calcium or magnesium in excess of 50 mg/L. Use of this set is in addition to digestion (using TNT890) when measuring total Cadmium.

Easy and safe handling

The innovative vial and reagent delivery makes tests easier to use than powder pillows or liquid reagents and significantly reduces spillage or contamination risk.

No reagent blank necessary

The high quality of TNTplus vials, tight reagent production controls, instrument calibration verification, and high instrument stability all combine to eliminate the need to run reagent blanks—saving you time and money!

Automatic method detection

The bar-coded vials allow for automatic method detection and measurement by the spectrophotometer, considerably increasing the speed of analysis. Errors are reduced as instrument averages 10 readings and rejects outliers from scratched, flawed, or dirty glassware.

Documented shelf life and COA

New 2D barcode details batch number and expiry date of reagents, which are documented along with the measurement result. An automatic warning is issued if expiry date has passed. Certificate of Analysis (COA) is available on RFID tag on the box.

Specifications

Instrument:	DR3900, DR6000, DR1900, DR2800, DR3800, DR5000
Parameter:	Cadmium
Platform:	TNTplus

Required Accessories

• Gloves, Disposable, Powder Free, Nitrile, Large (Item 2550503)