



Reagent powder pillows, DithiVer metals, pk/100

Product #: ZAR Price: Ships within 1 week **1261699** Contact Hach

Powder reagent in sealed package for the determination of Cadmium and Lead, Dithizone Method.

Specifications

| Method: | Dithizone |
|------------------|----------------|
| Number of tests: | 100 |
| Platform: | Powder Pillows |
| Quantity: | 100 /pk |

Required Accessories

• Gloves, Disposable, Powder Free, Nitrile, Large (Item 2550503)