



## Citrate Buffer Powder Pillows

<b>Product #:</b>	<b>2107669</b>
ZAR Price:	Contact Hach
Ships within 3 days	

Powder reagent in sealed packages for the determination of Manganese, cold periodate method. For 10 mL sample.

For determination of manganese by the periodate oxidation method. EPA approved for reporting wastewater analyses (digestion required). Hach method 8034. Sample size: 10 mL Pack of 100 Powder Pillows.

---

### Specifications

Method:	Periodate Oxidation
Number of tests:	100
Platform:	Powder Pillows
Quantity:	100 /pk