



## Sodium Periodate Powder Pillows, 10 mL, pk/100

**Product #:**

**2107769**

ZAR Price:  
Available

Contact Hach

 **Hazardous**

Items with this mark may be considered hazardous under some shipping conditions.

If necessary, we will change your selected shipping method to accomodate these items.

Powder reagent in sealed packages for the determination of Manganese, cold periodate method. For 10 mL sample.

For determination of manganese by the periodate oxidation method. EPA approved for reporting wastewater analyses (digestion required). Hach method 8034. Sample size: 10 mL Pack of 100 Powder Pillows.

---

### Specifications

Method:	Periodate Oxidation
Number of tests:	100
Platform:	Powder Pillows
Quantity:	100 /pk