



Sodium periodate powder pillows for manganese, 25 mL, pk/100

Product #:

98499

ZAR Price:

Contact Hach

Ships within 3 days

 **Hazardous**

Items with this mark may be considered hazardous under some shipping conditions.

If necessary, we will change your selected shipping method to accommodate these items.

Powder reagent in sealed package for the determination of Manganese, cold periodate method, and Iron titrations. For 25 ml sample.

For determination of manganese by the periodate oxidation method. EPA approved for reporting wastewater analyses (digestion required). Hach method 8034. Sample size: 25 mL Pack of 100 Powder Pillows.

Specifications

Method:	Periodate Oxidation
Number of tests:	100
Platform:	Powder Pillows
Quantity:	100 /pk