



Reagent set, manganese, HR, periodate oxidation, 25 mL



Product #:
ZAR Price:
Available

2243200
Contact Hach

 **Hazardous**

Items with this mark may be considered hazardous under some shipping conditions.

If necessary, we will change your selected shipping method to accommodate these items.

For determination of Manganese by the Periodate Oxidation Method, Hach Method 8034. EPA Approved for wastewater reporting (digestion required) - Federal Register, 44, (116) 34193 (June 14, 1979). Range: 0.1 - 20.0 mg/L Mn. Sample Size: 25 mL. Set Contains: Buffer Powder Pillows (Citrate type) and Sodium Periodate Powder Pillows for 100 tests.

Specifications

EPA compliant:	Yes
Instrument:	DR/3000, DR/2000
Method:	Periodate Oxidation
Number of tests:	100
Platform:	Powder Pillows
Range:	0.1 - 20.0 mg/L Mn

Required Accessories

- Gloves, Disposable, Powder Free, Nitrile, Small (Item 2550501)