



Manganese Reagent Set, Powder Pillows, periodate oxidation, 0.1-20.0 mg/L Mn, 25 mL sample

Product #:

2243200

 **Hazardous**

ZAR Price:

Contact Hach

Ships within 3 days

Items with this mark may be considered hazardous under some shipping conditions.

If necessary, we will change your selected shipping method to accomodate these items.

For determination of Manganese by the Periodate Oxidation Method, Hach Method 8034. EPA Approved for wastewater reporting (digestion required) - Federal Register, 44, (116) 34193 (June 14, 1979). Range: 0.1 - 20.0 mg/L Mn. Sample Size: 25 mL. Set Contains: Buffer Powder Pillows (Citrate type) and Sodium Periodate Powder Pillows for 100 tests.

Specifications

EPA compliant:	Yes
Instrument:	DR/3000, DR/2000
Method:	Periodate Oxidation
Number of tests:	100
Parameter:	Manganese
Platform:	Powder Pillows
Range:	0.1 - 20.0 mg/L Mn
Sample Volume:	25 mL