



## Manganese LR PAN (50)

Product #:

HPT291

ZAR Price:

Contact Hach

 **Hazardous**

Items with this mark may be considered hazardous under some shipping conditions.

If necessary, we will change your selected shipping method to accommodate these items.

For determination of Manganese by the PAN Method, using a 50-mm sample cell. Hach Method 10286. Ultra Low Range: 0.005 - 0.500 mg/L Mn. 50 tests/package.

**For determination of Manganese by the PAN method. Low Range 0.005 - 0.500 mg/L Mn.**