



Alkaline reagent powder pillows, HgEx reagent 3, pk/25

Product #:

2658448

Hazardous

ZAR Price: Ships within 3 days

Contact Hach

Items with this mark may be considered hazardous under some shipping conditions.

If necessary, we will change your selected shipping method to accomodate these items.

For mercury determination by the Cold Vapor Mercury Concentration Method. Hach Method 10065. Reagent is a component of reagent set 2658300.

Specifications

Method: Cold Vapor Concentration

Number of tests: 25

Platform: Powder Pillows

Required Accessories

• Gloves, Disposable, Powder Free, Nitrile, Large (Item 2550503)