



Complexing reagent powder pillows, HgEx reagent 6, pk/25

Product #:	2658748
ZAR Price:	Contact Hach
Ships within 2 weeks	

For mercury determination by the Cold Vapor Mercury Concentration Method. Hach Method 10065. Reagent is a component of reagent set 2658300.

Specifications

Method:	Cold Vapor Concentration
Number of tests:	25
Platform:	Powder Pillows

Required Accessories

- Gloves, Disposable, Powder Free, Nitrile, Large (Item 2550503)