



# MolyVer 6 molybdenum reagent AccuVac ampules, 25/pk

Product #:

2522025

🚹 Hazardous

**OBSOLETE ITEM** This item is no longer available.

Suggested replacements: • 2604100

hazardous under some shipping conditions. If necessary, we will change your selected shipping method to accomodate these items.

Items with this mark may be considered

For determination of HR Molybdate Molybdenum by the Mercaptoacetic Acid Method. Hach method 8036. Range: 0.2 - 40.0 mg/L Mo. Pack of 25 ampules. Also requires CDTA Solution, 0.4 M, Cat. No. 26154-36.

### Fast

AccuVac Ampules provide quick results. Just snap, fill and read - saving you valuable time.

#### Easy

No need to measure sample or reagents. The evacuated ampules draw in the correct sample volume, so anyone can obtain reliable results with minimal effort.

### Accurate

Optical quality evacuated vials assure precise measurements and reagent stability.

### Use them anywhere you need results.

Ampules are both reagent containers and photometric measurement cells. No need to transport reagents or samples.

## **Specifications**

Method:	Mercaptoacetic Acid
Number of tests:	25
Parameter:	Molybdenum
Platform:	Accuvac
Range:	0.2 - 40.0 mg/L Mo

## **Required Accessories**

• Gloves, Disposable, Powder Free, Nitrile, Large (Item 2550503)