



## Selenium Reagent Set (No air freight)

Product #: 2244200

 **Hazardous**

**OBSOLETE ITEM**  
This item is no longer available.

Items with this mark may be considered hazardous under some shipping conditions.

If necessary, we will change your selected shipping method to accommodate these items.

For determination of Selenium by the Diaminobenzidine Method. Hach Method 8194. Range: 0.01 - 1.00 mg/L Se. For air freight, order 2244201

---

### Specifications

Instrument:	All except DR 900, PC II
Method:	8194
Number of tests:	66 - 100
Platform:	Solution
Range:	0.01 - 1.00 mg/L Se

---

### Required Accessories

- Gloves, Disposable, Powder Free, Nitrile, Large (Item 2550503)