



## Reagent set, selenium (air freight configuration)

Product #: 2244201

 **Hazardous**

**OBSOLETE ITEM**  
This item is no longer available.

Items with this mark may be considered hazardous under some shipping conditions.

If necessary, we will change your selected shipping method to accomodate these items.

For determination of Selenium by the Diaminobenzidine Method. Hach Method 8194. Range: 0.01 - 1.00 mg/L Se. Reagents are packaged for air freight.

---

### Specifications

Method:	Diaminobenzidine
Number of tests:	66 - 100
Parameter:	Selenium
Range:	0.01 - 1.00 mg/L

---

### Required Accessories

- Gloves, Disposable, Powder Free, Nitrile, Large (Item 2550503)