



Dissolved oxygen 2 reagent powder pillows, 60 mL sample, pk/100

Product #: ZAR Price: Ships within 1 week 98299

Contact Hach



Items with this mark may be considered hazardous under some shipping conditions.

If necessary, we will change your selected shipping method to accomodate these items.

For determination of dissolved oxygen by the Winkler titration method. Pack of 100 powder pillows.

Specifications

Method: Platform: Winkler Powder Pillows

Required Accessories

• Gloves, Disposable, Powder Free, Nitrile, Large (Item 2550503)