



## Mercuric nitrate cartridge for digital titrator, 0.2256 N, 1 piece

**Product #:**

**1439301**

**ZAR Price:**

Contact Hach

**Ships within 3 days**



**Hazardous**

Items with this mark may be considered hazardous under some shipping conditions.

If necessary, we will change your selected shipping method to accommodate these items.

**Contains mercury. Dispose of according to local, state and federal regulations.**

Ready-to-use titrant solution for chloride determination by Digital Titrator titration with mercuric nitrate, using diphenylcarbazone indicator. Hach method 8206. Approximately 40-100 tests.

### Accurate results - when and where you need them

Digital Titrator methods use a high-precision dispensation device and cartridges of concentrated titrant to allow analyses in the lab, the plant or the field.

### Fast and Flexible

Interchangeable cartridges allow you to change titrants in seconds. No cleaning or rinsing burets.

### Versatile

More than 50 parameter/range combinations are available

### Economical

Concentrated titrants allow 40-100 typical titrations per cartridge.

---

## Specifications

Method:	Mercuric Nitrate
Number of tests:	40 - 100
Parameter:	Chloride
Platform:	Digital Titrator
Range:	10 - 160 mg/L Cl <sup>-</sup>
Titant:	0.2256 N Hg(NO <sub>3</sub> ) <sub>2</sub>

---

## Required Accessories

- Gloves, Disposable, Powder Free, Nitrile, Large (Item 2550503)