



Mercuric nitrate cartridge for digital titrator, 2.256N

Product #:

ZAR Price: Ships within 3 days 92101

Contact Hach

А н

Hazardous

Items with this mark may be considered hazardous under some shipping conditions.

If necessary, we will change your selected shipping method to accomodate these items.

Contains mercury. Dispose of according to local, state and federal regulations.

For titration of chloride, 100 - 8000 mg/L Cl

Specifications

Method: 8206

Number of tests: 40 - 100
Parameter: Chloride

Platform: Digital Titrator

Range: 100 - 8000 mg/L Cl⁻
Titrant: 2.256 N Hg(NO₃)₂

Required Accessories

• Gloves, Disposable, Powder Free, Nitrile, Large (Item 2550503)