



Digital titrator cartridge, mercuric nitrate, 2.570 N

Product #: 2393701
ZAR Price: Contact Hach
Ships within 3 days

 **Hazardous**

Items with this mark may be considered hazardous under some shipping conditions.

If necessary, we will change your selected shipping method to accommodate these items.

Contains mercury. Dispose of according to local, state and federal regulations.

For determination of Salinity by Digital Titrator titration. 2.570 N. Digits x 0.1 = ppt Salinity in 2 mL sample.

Accurate results - when and where you need them

Digital Titrator methods use a high-precision dispensation device and cartridges of concentrated titrant to allow analyses in the lab, the plant or the field.

Fast and Flexible

Interchangeable cartridges allow you to change titrants in seconds. No cleaning or rinsing burets.

Versatile

More than 50 parameter/range combinations are available

Economical

Concentrated titrants allow 40-100 typical titrations per cartridge.

Specifications

Instrument:	Digital Titrator
Method:	Mercuric Nitrate
Number of tests:	40 - 100
Parameter:	Sulfite
Platform:	Digital Titrator
Range:	0 - 100 ppt Salinity
Titrant:	2.570 N Hg(NO ₃) ₂

Required Accessories

- Gloves, Disposable, Powder Free, Nitrile, Large (Item 2550503)