



Phenylarsine Oxide (PAO) Digital Titrator Cartridge, 0.0451 N

Product #:	2142001
ZAR Price:	Contact Hach
Available	

For chlorine determination by Digital Titrator titration. 0.00451 N. 1 digit = 1 mg/L Cl_2 in 200 mL sample, Digits x 0.00000564 = meq.

Accurate results - when and where you need them

Digital Titrator methods use a high-precision dispensation device and cartridges of concentrated titrant to allow analyses in the lab, the plant or the field.

Fast and Flexible

Interchangeable cartridges allow you to change titrants in seconds. No cleaning or rinsing burets.

Versatile

More than 50 parameter/range combinations are available

Economical

Concentrated titrants allow 40-100 typical titrations per cartridge.

Specifications

Method:	Amperometric Titration
Platform:	Digital Titrator

Required Accessories

- Gloves, Disposable, Powder Free, Nitrile, Large (Item 2550503)