



## Plan Objective 10X for Leica DM500

Product #:

2945028

 **Hazardous**

**OBSOLETE ITEM**

This item is no longer available.

Items with this mark may be considered hazardous under some shipping conditions.

If necessary, we will change your selected shipping method to accommodate these items.

Each objective is labeled with its magnification, for instance 5x or 100x. However, the magnification of the objective alone does not determine the overall magnification of the microscope. This results from the objective magnification multiplied by the eyepiece magnification (for tube lens 1x). Example: 40x objective x 10x eyepieces = 400x overall magnification. It should be noted, however, that the higher the objective magnification, the lower the visible object field.