



## M COLIBLUE24 Broth ampules (pk/50)

<b>Product #:</b>	<b>2608450</b>
ZAR Price:	Contact Hach
Ships within 3 days	

In 24 hours, Hach's m-ColiBlue24 Broth give you simultaneous results for total coliforms and *E. coli*. Hach's ready-to-use media eliminates measuring, mixing, and autoclaving steps necessary to prepare dehydrated media. Glass ampules (26084-20), plastic ampules (260845-00), 100mL bottles (26084-42) of m-ColiBlue24 and agar plates (28052-15), afford maximum shelf life and ease of use for pre-measured and prepared media. Enumerate total coliforms and *E. coli* on one petri dish. Differentiate colonies easily - red and blue indicate total coliforms and blue specifies *E. coli*. Get optimal recovery of stressed and injured organisms.

Use m-ColiBlue24 Broth to accurately monitor and evaluate total coliform and *E. coli* in drinking water; wastewater; bottled water; beverages; surface, ground, and well water; and ultrapure, chemical processing and pharmaceutical processing waters.

m-ColiBlue24 should only be used with pads manufactured by the following companies/brands: Pall Corporation (e.g., Sterile Petri Dishes with Absorbant Pads), Private Label/OEM products manufactured by Pall Gelman, Sartorius (cellulose and glass fiber pads), Sartorius BioSart products with the integrated Petri dish, and Millipore EO sterilized pads. No other pad is compatible with the product, and the use of any pad other than those recommended may result in unreliable, inaccurate data.

With ready-to-use m-ColiBlue24 Broth, you can:

**Simultaneously detect total coliforms and *E. coli* - in just 24 hours**

**Simple handling with pre-portioned selective medium**

**Cost-effective - no fluorescent lamps or colour comparators required**

**Superior sensitivity down to 1CFU/100 mL**

---

### Specifications

Material:	Plastic
Media Type:	m-Colibblue24 Broth
Quantity:	50 /pk