



Tubes, Micro Dist Cyanide-1, NaOH trapping solution, pk/21

Product #:

A17001

ZAR Price:

Contact Hach

Ships within 3 days

 **Hazardous**

Items with this mark may be considered hazardous under some shipping conditions.

If necessary, we will change your selected shipping method to accommodate these items.

Tubes for Cyanide using the Micro Dist Distillation System. Cyanide method is USEPA Approved for NPDES and NPDWR

Save time and money with disposable Micro Dist tubes that contain the NaOH trapping solution!

Increase reliability with press-fit joints that will never leak

Eliminates bulky, expensive glassware

Accommodates liquid and solid samples

Less Hazardous for users!