



## **CDC749 Conductivity cell, 2-pole, platinised, microsample, d=4 mm, MAB6 plug (Radiometer Analytical)**

<b>Product #:</b>	<b>E61M014</b>
ZAR Price:	Contact Hach
Ships within 2 weeks	

For use in glass tubes.

### **The right conductivity cell for your application.**

It is important to choose a conductivity cell with the right construction and geometry for your particular application and working conditions. The Radiometer Analytical line offers a variety of conductivity cells to meet a large number of applications.

---

### **Specifications**

Compatible Instruments:	MeterLab CDM210, CDM230
Connector:	MAB6
Parameter:	Conductivity
Type:	Lab, 2 pole Pt
Warranty:	12 months