



Chloride Standard Solution, 1000 mg/L, 500 mL

Product #: 18349
ZAR Price: Contact Hach
Ships within 3 days

Ion Selective Electrode (ISE) standard solutions for laboratory analyses

ISE Standard Solution. Chloride, 1000 mg/L ± 10 mg/L Cl^- (NIST). 500 mL bottle

Traceable to NIST standard reference materials

Always performing measurement results you can trust is essential. The quality of the standard solutions that are used are key.

Design for use with Hach meters portfolio and other meters.

The basis for accurate ISE measurements is good care and calibration of the Ion Selective Electrode.

Reliable packaging

Supplied in an airtight packaging ensuring guaranteed shelf-life.

Specifications

Concentration:	1000 mg/L ± 10 mg/L Cl^-
Concentration 2:	0.02820 mol/L ± 0.00028 mol/L
Description:	Chloride Standard Solution as Cl^- (NIST)
Package Type:	Bottle
Parameter:	Chloride
Platform:	Solution
Quantity:	500 mL
Special Remarks:	Observed shelf life after opening, following recommended storage and use: 3 months
Storage conditions:	10 - 25 °C (protect from light)

Required Accessories

- Chloride Ionic Strength Adjustor (ISA) Powder Pillows, pk/100 (Item 2318069)