



## Ammonia TNTplus Vial Test, HR (2-47 mg/L $\text{NH}_3\text{-N}$ ), 25 Tests

**Product #:**  
**ZAR Price:**  
Ships within 1 week

**TNT832**  
Contact Hach



Items with this mark may be considered hazardous under some shipping conditions.

If necessary, we will change your selected shipping method to accomodate these items.

### Expert High Range Ammonia test made simple.

Ammonia TNTplus chemistry and your Hach spectrophotometer are engineered to simplify water analysis for accurate results, everytime.

For determination of Ammonia Nitrogen by the Salicylate Method. EPA compliant.

#### Easy and safe handling

The innovative vial and reagent delivery makes tests easier to use than powder pillows or liquid reagents and significantly reduces spillage or contamination risk.

#### No reagent blank necessary

The high quality of TNTplus vials, tight reagent production controls, instrument calibration verification, and high instrument stability all combine to eliminate the need to run reagent blanks—saving you time and money!

#### Automatic method detection

The bar-coded vials allow for automatic method detection and measurement by the spectrophotometer, considerably increasing the speed of analysis. Errors are reduced as instrument averages 10 readings and rejects outliers from scratched, flawed, or dirty glassware.

#### Documented shelf life and COA

New 2D barcode details batch number and expiry date of reagents, which are documented along with the measurement result. An automatic warning is issued if expiry date has passed. Certificate of Analysis (COA) is available on RFID tag on the box.

---

### Specifications

EPA compliant:	Yes
Instrument:	DR3900, DR6000, DR1900, DR2800, DR3800, DR5000
Method:	10205
Number of tests:	25
Parameter:	Ammonia, Nitrogen
Platform:	TNTplus
Range:	2 - 47 mg/L $\text{NH}_3\text{-N}$
Storage conditions:	2 - 8 °C (keep refrigerated)
Truecal:	Yes

---

## Required Accessories

- Gloves, Disposable, Powder Free, Nitrile, Large (Item 2550503)