



Volatile Acids TNTplus Vial Test (50-2,500 mg/L), 25 Tests

Product #:

TNT872

ZAR Price:

Contact Hach

Ships within 3-5 weeks



Hazardous

Items with this mark may be considered hazardous under some shipping conditions.

If necessary, we will change your selected shipping method to accommodate these items.

Expert Volatile Acids test made simple.

Volatile Acids TNTplus chemistry and your Hach spectrophotometer are engineered to simplify water analysis for accurate results, everytime.

For determination of Volatile Acids by the Esterification method.

Simple test allows quick process control - No need for titrations or involved water bath digestion.

Easy and safe handling

The innovative vial and reagent delivery makes tests easier to use than powder pillows or liquid reagents and significantly reduces spillage or contamination risk.

No reagent blank necessary

The high quality of TNTplus vials, tight reagent production controls, instrument calibration verification, and high instrument stability all combine to eliminate the need to run reagent blanks—saving you time and money!

Automatic method detection

The bar-coded vials allow for automatic method detection and measurement by the spectrophotometer, considerably increasing the speed of analysis. Errors are reduced as instrument averages 10 readings and rejects outliers from scratched, flawed, or dirty glassware.

Documented shelf life and COA

New 2D barcode details batch number and expiry date of reagents, which are documented along with the measurement result. An automatic warning is issued if expiry date has passed. Certificate of Analysis (COA) is available on RFID tag on the box.

Specifications

Digestion Required:	Yes
EPA compliant:	Yes
Instrument:	DR3900, DR6000, DR1900, DR2800, DR3800, DR5000
Method:	10240
Number of tests:	25
Parameter:	Volatile Acids
Platform:	TNTplus
Range:	50 - 2,500 mg/L Acetic Acid

Required Accessories

- Gloves, Disposable, Powder Free, Nitrile, Large (Item 2550503)