



BOD Seed Inoculum, Polyseed, pk/50 capsules

Product #: 2918700
ZAR Price: Contact Hach
Ships within 3 days

BOD Seed Inoculum for use in Biochemical Oxygen Demand (BOD5) test procedures. A blend of specialized microbial cultures in an easy-to-use capsule. Designed to provide a uniform standard for the degradation of both industrial and municipal waste in BOD5 analyses. Package of 50 capsules. International Laboratory Supply (InterLab), LTD

- * Consistent seed source: same bacteria, same type of bacteria, same amount of bacteria blend
- * Known bacterial content is free from pathogens
- * Save and convenient pre-measured dry product in easy to use capsules
- * Two years shelf life

Specifications

| | |
|---------------|----------------------------|
| Application: | Dilution water preparation |
| Description: | Seed Source |
| Formula: | Seed Source |
| Method: | 8043 |
| Package Type: | Capsules |
| Parameter: | Biochemical Oxygen Demand |
| Quantity: | 50 pcs |

Required Accessories

- Gloves, Disposable, Powder Free, Nitrile, Large (Item 2550503)