



## VDK (vicinal diketones) TNTplus Vial Test (0.015-0.5 mg/kg)

**Product #:**

**TNT819**

**ZAR Price:**

Contact Hach

No availability shown

 **Hazardous**

Items with this mark may be considered hazardous under some shipping conditions.

If necessary, we will change your selected shipping method to accommodate these items.

**Avoid off-flavors in your beer and safeguard your product quality by measuring VDK (vicinal diketones).**

For use with your Hach DR6000 spectrophotometer. TNT819 offers you the ability to test for the presence of VDK with less hassle and chemical mess.

Additionally, these vials provide high quality results free from potential external influencers. They are designed to simplify your brewing analysis, with accurate results, every time.

Features include:

**Ready-to-use, high-precision pre-dosed cuvette tests**

**Reduction of errors due to bar-coded vials**

**RFID technology to document the shelf life of the cuvettes and ensure traceability of your analyses**

**Helpful Package Design**

---

### Specifications

Instrument:	DR6000, DR3900, DR5000, DR3800
Method:	10276
Number of tests:	25
Parameter:	Vicinal diketones (VDK)
Platform:	TNTplus
Range:	0.015 - 0.5 mg/kg Diacetyl
Storage conditions:	2 - 8 °C

---

### Required Accessories

- Gloves, Disposable, Powder Free, Nitrile, Large (Item 2550503)