



**Be Right™**



## **Flask, Erlenmeyer, Polymethylpentene Capacity 500 mL, 4/ pk**

**Product #:**

**2089879**

**ZAR Price:**

Contact Hach

**Ships within 1 week**

Polymethylpentene flask with polypropylene screwcap features clarity, chemical resistance, strength, and autoclavability (Do not engage threads during cycle). Useful for titrations, preparation and storage of media.

---

### **Specifications**

Capacity:	500 mL
Material:	Polymethylpentene
Type:	Erlenmeyer
Units:	4/pk