



Bottle, wash, diamond alert, 500 mL, Variety 6/pk

Product #:	2613049
ZAR Price:	Contact Hach
Ships within 1 week	

Wide-mouthed and color-coded wash bottles provide safe and positive identification of common laboratory solutions. Standard Diamond Alert symbols indicate type of hazard contents may represent. Polyethylene bottles with color-coded polypropylene screw caps and spray-tip assembly. Variety 6-pack includes one bottle each for water, acetone, methanol, isopropanol, toluene, and ethanol.

Specifications

Capacity:	500 mL
Color:	Transparent
Material:	LDPE
Quantity:	6 pk
Stem:	Blue Cap/Stem
Unit:	6 pk