



Beaker, Griffin, Low Form, Glass, 250 mL, 12/pk

| | |
|------------|--------------|
| Product #: | 50076H |
| ZAR Price: | Contact Hach |
| Available | |

Griffin style (low form) beakers are the industry standard for a wide variety of lab applications. Graduations are high quality screen printing. Bull nose lip and molded pour spout minimize spillage in liquid transfer. Autoclavable, non-toxic and biologically inert. Volumes and graduations are approximate.

Specifications

| | |
|----------------|--------------------------|
| Body material: | Glass |
| Capacity: | 250 mL |
| Material: | Glass, Griffin, Low Form |
| Units: | 12 /pk |