



Replacement: TUBING, TYGON(R), 3/8" ID X 5/8" OD

| | |
|-------------------|----------------|
| Product #: | 2943221 |
| ZAR Price: | Contact Hach |
| Available | |

Replacement Tygon Tubing, 3/8" interior diameter X 5/8" outside diameter for Xenosep SPE Systems Apparatus.

Specifications

| | |
|-----------------|--------|
| Body material: | Tygon |
| Diameter: | 5/8 in |
| I. D.: | 3/8 in |
| Length: | 3 ft |
| Wall Thickness: | 1/4 in |