



Test tube, rack, polyethylene, 25 mm (O.D.), 24 holes

Product #:	2497903
ZAR Price:	Contact Hach
Available	

Twenty to ninety-place rack features a special blend of polypropylene which sinks in water, assuring non-floating racks when immersed. Autoclavable. Stackable and autoclavable. Overall dimensions: 24.6 x 10.7 x 6.4 cm (9.7 x 4.2 x 2.5").

Specifications

Diameter:	25 mm
Holes:	24
Quantity:	1 each
Test Tube Arrangement:	3 x 8