



Digital thermometer, NIST

Product #:	2630600
ZAR Price:	Contact Hach
Ships within 1 week	

Light weight and inexpensive, this K type thermometer provides extremely accurate and reliable temperature measurement. Features relative minimum and maximum temperature, hold, auto power off, and touch tone buttons. This is the meter only. For use with probe 2630601 (sold separately).

Specifications

Length:	133 mm
Quantity:	1 each
Temperature Range:	-328 °F - 2498 °F
Thermometer Divisions:	0.1 °F

Required Accessories

- Stainless Steel Replacement K-Type Probe (Item 2630601)