



Thermometer, Non-Mercury, 0 - 220 °F

| | |
|-------------------|----------------|
| Product #: | 2676200 |
| ZAR Price: | Contact Hach |
| Available | |

Compact and durable, non-mercury thermometers enclosed in a protective nylon case. These thermometers are designed for ruggedness and come in many standard ranges. Suitable for all field and portable uses where an unprotected thermometer may break.

Specifications

| | |
|-------------------------|---------------|
| Quantity: | 1 each |
| Temperature Range: | 0 °F - 220 °F |
| Thermometer Divisions: | 2 °F |
| Thermometer Fluid Type: | Non-Mercury |