



Thermometer, Enviro-Safe, Non-Mercury, 0 to 300°F

Product #:	2635704
ZAR Price:	Contact Hach
Available	

Enviro-Safe non-Mercury thermometers feature an exclusive biodegradable liquid and dye. Easy to read, aesthetically distinctive due to dark green liquid and permanently fused markings. Each thermometer, individually calibrated against NIST traceable equipment, identified by unique serial number for ISO

9000 or other quality systems.

Specifications

Length:	305 mm
Temperature Range:	0 °F - 300 °F
Thermometer Divisions:	2 °F