



## Test kit, surface water

**Product #:**

**2559800**

**ZAR Price:**

**Contact Hach**

**Call for ship date**

 **Hazardous**

Items with this mark may be considered hazardous under some shipping conditions.

If necessary, we will change your selected shipping method to accomodate these items.

The ideal tool for source water monitoring and wastewater treatment plants, this kit uses drop count titration and color disc methods to provide economical, easy-to-use tests for ammonia, chlorine, pH, nitrate, dissolved oxygen, phosphorus, and temperature. The rugged, lightweight kit includes the convenient pH Pocket Pal Tester and all necessary apparatus and reagents for testing.

**Ideal tool for both industrial and municipal water treatment plants**

**Uses drop count titration and color disc methods to provide economical, easy-to-use tests for ammonia, chlorine, pH, titrate, dissolved oxygen, phosphorus, and temperature**

**Rugged, lightweight kit includes the convenient pH Pocket Pal™ Tester and all necessary apparatus and reagents**

---

### Specifications

Application:	Surface Water
Number of tests:	100
Platform:	Color Disc