



Free & Total Chlorine Test Kit (Model CN-70)

Product #: 1454200
ZAR Price: Contact Hach
Ships within 3 days

Simple, low-cost color disc method is the most accurate visual test kit method

Simple, low-cost color disc method is the most accurate visual test kit method

Achieve fast results with an easy-to-carry test kit

Specifications

Case Style:	D
Footnote:	*mg/L unless otherwise noted, ppb = $\mu\text{g/L}$, ppm = mg/L.; gpg = grains per gallon; 1 gpg = 17.1 mg/L or 17.1 ppm.
Measuring range:	Dual Range
Method:	Color disc/DPD
Number of tests:	100 Free, 100 Total
Parameter:	Chlorine (Free and Total)
Platform:	Color Disc
Range:	0.02 - 0.68 mg/L Cl_2
Range 2:	0 - 3.4 mg/L Cl_2
Smallest increment:	0.02
	0.1
Test Kit Model:	CN-70 ²