



## Nitrogen, Nitrite Colour Disc Test Kit, Model NI-6

**Product #:** 224000  
**ZAR Price:** Contact Hach  
**Ships within 3 days**

Utilising the simple colour disc method/ferrous sulfate, this kit measures Nitrogen, Nitrite as  $\text{NO}_2^-$  and contains 100 tests. Range (mg/L): 0-8/0-1600.

**Simple, low-cost colour disc method is the most accurate visual test kit method**

**Kit contains all required reagents and apparatus in a rugged carrying case for analysis in the field**

---

### Specifications

Case Style:	D
Measuring range:	Single Range
Method:	Colour disc/  Ferrous sulphate
Number of tests:	100
Parameter:	Nitrogen, Nitrite
Platform:	Color Disc
Range:	0 - 8/0-1600 mg/L $\text{NO}_2^-$
Smallest increment:	2  40
Test Kit Model:	NI-6
What's included?:	Includes colour disc, comparator box and tubes, reagents, deionised dilution water, instructions and carrying case.

---

### What's included?

Includes colour disc, comparator box and tubes, reagents, deionised dilution water, instructions and carrying case.