



## Hardness (Total & Calcium) Test Kit, Model HA-4P

**Product #:**

**145700**

**ZAR Price:**

**Contact Hach**

**Ships within 3 days**

 **Hazardous**

Items with this mark may be considered hazardous under some shipping conditions.

If necessary, we will change your selected shipping method to accommodate these items.

Hardness, Total & Calcium drop count titration/EDTA test kit, range 17.1 - 342 mg/L  $\text{CaCO}_3$ , 100 tests.

The kit contains all required reagents and apparatus in a rugged carrying case for analysis in the field.

**Utilises drop count titration - a simple, economical titration method**

---

### Specifications

Case Style:	D
Measuring range:	17.1 - 342 mg/L $\text{CaCO}_3$
Method:	Drop Count Titration/EDTA
Number of tests:	100
Parameter:	Hardness, Total & Calcium - As $\text{CaCO}_3$
Platform:	Drop Count Titration
Smallest increment:	17.1 mg/L
Test Kit Model:	HA-4P
What's included?:	Includes five reagents, bottle, measuring tube, instruction sheet, and carrying case.

---

### What's included?

Includes five reagents, bottle, measuring tube, instruction sheet, and carrying case.