



Chlorine (total) test kit, model CN-21P, 10 - 200 mg/L, 100 tests

Product #:

2444400

ZAR Price:

Contact Hach

Ships within 3 days



Hazardous

Items with this mark may be considered hazardous under some shipping conditions.

If necessary, we will change your selected shipping method to accomodate these items.

Method / Chemistry: Drop count titration / Thiosulfate

Chlorine Test Kits are manufactured by Hach Company to ensure compliance with stringent quality control standards. Using the Drop Count Titration method, the CN-21P Chlorine Test Kit is used to determine high concentrations of chlorine in sanitizing applications. 10 - 200 mg/L; contains 100 tests.

Simple, economical titration method

Kit contains all required reagents, solutions, and glassware

Specifications

Case Style:	B
Footnote:	*mg/L unless otherwise noted, ppb = $\mu\text{g/L}$, ppm = mg/L.; gpg = grains per gallon; 1 gpg = 17.1 mg/L or 17.1 ppm.
Measuring range:	10 - 200 mg/L
Method:	Drop count titration/ Thiosulfate
Number of tests:	100
Parameter:	Chlorine, mid and high range - As total Cl_2
Platform:	Drop Count Titration
Smallest increment:	10
Test Kit Model:	CN-21P