



Test kit, alkalinity, model AL-DT

Product #:

2063700

ZAR Price:

Contact Hach

Ships within 3 days



Hazardous

Items with this mark may be considered hazardous under some shipping conditions.

If necessary, we will change your selected shipping method to accommodate these items.

Caused by the presence of carbonates, bicarbonates, hydroxides, and other dissolved salts, alkalinity is important to food and beverage processors, boiler operators, and aquatic biologists. Determine phenolphthalein and total alkalinity by visual endpoint titration. Using the Digital Titrator, this contains 100 tests with a range of 10-4,000 mg/L.

Digital Titrator and standardized titration cartridges provide laboratory accuracy for field titrations

Kit contains all required reagents and apparatus in a rugged carrying case for analysis in the field

Specifications

Case Style:	F
Footnote:	*mg/L unless otherwise noted, ppb = $\mu\text{g/L}$, ppm = mg/L.; gpg = grains per gallon; 1 gpg = 17.1 mg/L or 17.1 ppm. 1EPA-accepted for wastewater. Check with your regulatory agency.
Method:	Digital Titrator/Sulfuric acid
Number of tests:	100
Parameter:	Alkalinity - (P) & Total (MO) - As CaCO_3
Platform:	Digital Titrator
Range:	10 - 4,000 mg/L CaCO_3
Smallest Increments Steps:	0.1-10
Test Kit Model:	AL-DT