



## Test kit, chloride, model CDS-DT

**Product #:**

**2580600**

**ZAR Price:**

Contact Hach

**Ships within 1 week**

 **Hazardous**

Items with this mark may be considered hazardous under some shipping conditions.

If necessary, we will change your selected shipping method to accommodate these items.

The single-parameter Model CDS-DT uses the Hach Digital Titrator with Silver Nitrate Titrant Cartridges to determine chloride concentrations. Kit includes Digital Titrator, reagents and labware, manual and carrying case. Range: 10-10,000 mg/L. Contains approximately 100 tests.

**Hach Digital Titrator offers higher levels of precision and accuracy than drop count procedures, with reproducibility comparable to a buret titration**

**Rugged titration device is ideal for field use, with no glass burets required**

---

### Specifications

Case Style:	F
Method:	Digital Titrator/  Silver nitrate
Number of tests:	100
Parameter:	Chloride - As $\text{Cl}^-$ , 10-10,000 mg/L $\text{Cl}^-$
Platform:	Digital Titrator
Range:	10 - 10,000 mg/L $\text{Cl}^-$
Range 2:	to
Smallest increment:	10
Test Kit Model:	CDS-DT

---

### Required Accessories

- Pipette, variable, volume 1.0-5.0 mL (Item BBP065)
- Gloves, Disposable, Powder Free, Nitrile, Medium (Item 2550502)
- Gloves, Disposable, Powder Free, Nitrile, X-Large (Item 2550504)
- Gloves, Disposable, Powder Free, Nitrile, Small (Item 2550501)
- Gloves, Disposable, Powder Free, Nitrile, Large (Item 2550503)